




Solutions Guide for Sanger Sequencing Poor Quality

Description

Poor quality results are typically due to an inefficient reaction. In general, some or all of the trace before position 500 is close to background signal, which makes base calling in these regions unreliable. Furthermore, results are scored as “poor quality” if they fail to pass quality control and cannot be well classified as another failure type. Please note that the distinction between “poor quality” and “no priming” results is sometimes arbitrary; however, “poor quality” data may provide some useful data whereas “no priming” usually does not.



*The average signal intensity of each base (G, A, T, C) is displayed at the bottom of the online trace viewer.

Troubleshooting

There are several possible causes for a poor quality result. Explanations and solutions are provided in the table below. The list of potential causes is not exhaustive but represents the most common issues. For more assistance, please contact Technical Support team at +49-341-520122-41 or sanger.europe@genewiz.com.

POSSIBLE CAUSE	SOLUTION
<p>Primer concentration is not in the range of our guidelines. The primer concentration could be too low or high.</p>	<ul style="list-style-type: none"> Check that the primer concentration is in the correct range (If not using a GENEWIZ universal primer) <ul style="list-style-type: none"> Premix reaction: Add 5 µL of 5 pmol/µL primer to your DNA samples. Pre-Defined reaction: Send your primer separately at 10 pmol/µL (20 µL for tube reactions, 60 µL for reactions submitted in plates).



Solutions Guide for Sanger Sequencing Poor Quality


POSSIBLE CAUSE	SOLUTION
<p>Concentration of DNA template is not in the range of our guidelines. The DNA template concentration could be too low or high.</p>	<ul style="list-style-type: none"> • Double check that the DNA concentration is in the correct range. Review our Sample Submission Guidelines. • For PCR purified products please apply the appropriate method to determine concentration: <ul style="list-style-type: none"> ○ For PCR products purified with a DNA-binding matrix (e.g. column or beads), a spectrophotometer can be used to measure the concentration; however, accuracy declines below 10 ng/μL. For best results, measure the concentration of the eluate before diluting it. ○ For PCR products purified with an enzymatic cleanup protocol (e.g. ExoSAP-IT), a spectrophotometer cannot be used for accurate DNA quantitation. Reaction components, such as degraded primers and nucleotides, will absorb UV light and inflate the calculated DNA concentration. Instead, use the band intensity on an agarose gel relative to that of mass standards to estimate DNA concentration. Alternatively, a fluorometer, which employs a dye that specifically binds dsDNA, can provide accurate quantitation. • For unpurified PCR products the following can be performed to verify the concentration: <ul style="list-style-type: none"> ○ Run your product on an agarose gel to confirm that it produces a strong single band of expected size. If not, optimize your PCR conditions. • For additional tips on how to optimize PCR for sequencing, see our Sanger Tips & Tricks series. • We can repeat reactions with an alternative protocol using a less or more template volume, this could improve the data quality (applies only to our predefined service type).
<p>Plasmid is greater than 20 kb. Larger plasmids can be difficult to be sequenced efficiently and rolling circle amplification is less efficient on larger templates.</p>	<ul style="list-style-type: none"> • For plasmids greater than 20 kb, ensure you submit highly pure samples at 100-150 ng/μL. • Use PCR to amplify the region of interest and submit PCR products for sequencing.



GENEWIZ®
Solid science. Superior service.



Solutions Guide for Sanger Sequencing Poor Quality

POSSIBLE CAUSE	SOLUTION
<p>Secondary structures. Template forms secondary structure</p>  <p>Template</p>	<ul style="list-style-type: none"> • Try sequencing with a primer that is further upstream or downstream of the current primer location or try sequencing from the other direction. Sometimes changing the primer location can help sequence through a difficult region. • We can repeat the reaction using GENEWIZ's "Alternative Protocol", designed for difficult templates. This protocol helps to denature secondary structure that impedes polymerization.
<p>Low copy plasmid used. In case Bacterial Colonies or Glycerol Stocks were send these plasmids may not provide sufficient input DNA for efficient amplification.</p>	<ul style="list-style-type: none"> • Use a vector with a high copy origin of replication. • Prep the DNA and submit purified plasmids. • Use PCR to amplify the region of interest and submit PCR products for sequencing.
<p>EndA+ strain used. In case Bacterial Colonies or Glycerol Stocks were send, please note that certain strains of <i>E. coli</i> (e.g. BL21, Stbl3) contain a non-specific endonuclease in the periplasmic space that can cleave plasmid DNA when cells are lysed, leading to inefficient amplification.</p>	<ul style="list-style-type: none"> • Use a strain with the <i>endA</i> mutation (e.g. DH5α, TOP10). • Prep the DNA and submit purified plasmids. • Use PCR to amplify the region of interest and submit PCR products for sequencing.
<p>Rich media used. When sending Glycerol stocks or liquid bacterial cultures, the media used could interfere with the sequencing reaction. Certain media (e.g. TB, SOB, 2YT) contain high salt concentration or ingredients that can inhibit polymerase activity during rolling circle amplification.</p>	<ul style="list-style-type: none"> • Use LB medium. • Submit Bacterial Colonies.