

MagicPure® Viral DNA/RNA Kit

Cat. No. EC301

Storage: Proteinase K at -20°C for two years; Magnetic Virus Beads at 2-8°C for one year (protect from freezing); others at room temperature (15°C-25°C) for one year.

Description

MagicPure® Viral DNA/RNA Kit utilizes a unique lysis buffer to lyse virus and release DNA / RNA. The released DNA / RNA is effectively purified after specifically binding to silica magnetic beads. It is suitable for isolating viral DNA/RNA from up to 200 µl of plasma, serum, whole blood, tissue homogenate, cell-free body fluid, nasopharyngeal or oropharyngeal aspirate/ wash, bronchoalveolar lavage fluid (BALF), tracheal aspirate, sputum, nasopharyngeal or oropharyngeal swab and animal cell culture supernatant. The isolated DNA/RNA with high purity can be applied in PCR, RT-PCR, qPCR, and qRT-PCR, etc. This kit is compatible with high-throughput magnetic-rod nucleic acid extractor.

Highlights

- Simple and fast, no centrifugation required
- High yield, high purity

Kit Contents

Component	EC301-01/11 (50 rxns)
Binding Buffer 19(BB19)	18 ml
Clean Buffer 19(CB19)	18 ml
Wash Buffer 19(WB19)	12 ml
RNase-free Water	10 ml
Carrier RNA(1 µg/µl)	10 ml
Proteinase K (20 mg/ml)	1 ml
Magnetic Virus Beads	1 ml
Magnetic Stand(16 hole)	1 each/-

Sample requirements

- Store at 4°C for no more than 72 hours; at -70°C for long term storage
- Avoid repeated freezing and thawing
- Swab samples should only be collected with synthetic tip swabs (such as polyester or Dacron®) with aluminum or plastic shafts.

Procedures

Before starting, add 6 ml isopropanol to BB19, and add 12 ml and 48 ml anhydrous ethanol to CB19 and WB19, respectively.

1. Sample processing

• Liquid samples

(a) Add 20 µl Proteinase K and 300 µl BB19 to a sterile 1.5 ml microcentrifuge tube.

For multiple samples, mix Proteinase K and BB19 by the above ratio and separate the mixture into aliquots of 320 µl.

(b) Add 200 µl of sample to the microcentrifuge tube (If the sample volume is less than 200 µl, adjust the sample volume to 200 µl with 1×PBS or 0.9% NaCl.). Mix by vortexing for 5 seconds.

(c) Add 20 µl of Magnetic Virus Beads (vortex to ensure the beads are kept in suspension before use) and mix by vortexing for 30 seconds. Incubate at room temperature for 10 minutes and invert the tube 3-5 times during the incubation.



- **Solid samples (e.g., swabs):**

- (a) Place a swab and the entire storage buffer into a sterile 1.5 ml microcentrifuge tube and cut off the swab tip.
- (b) Add 450 μ l BB19 and 20 μ l Proteinase K. Mix by vortexing.
- (c) Incubate at 56°C for 20 minutes and vortex 3-5 times during the incubation.
- (d) Remove the swab tip, centrifuge briefly, and add 20 μ l of Magnetic Virus Beads (vortex to ensure the beads are kept in suspension before use) and mix by vortexing for 30 seconds. Incubate at room temperature for 10 minutes and invert the tube 3-5 times during the incubation.

- **For viscous liquids such as sputum, refer to " Solid Samples".**

2. Place the microcentrifuge tube on the magnetic stand until the beads are pelleted against the magnet. Discard the supernatant carefully (avoid pipetting the beads).

Recommendations on magnetic beads pelleting: After placing the microcentrifuge tube on the magnetic stand, gently turn the tube left and right until the beads are attached to the magnet. Gently invert the magnetic stand 2-3 times to ensure the beads on tube cap attached to the magnet. Incubate for 1 minute.

3. Remove the microcentrifuge tube from the magnetic stand and add 500 μ l of CB19 (make sure anhydrous ethanol has been added) to the tube. Mix by vortexing for 15 seconds and proceed to magnetic beads pelleting. Discard the supernatant carefully and avoid pipetting the beads.

4. Remove the microcentrifuge tube from the magnetic stand and add 500 μ l of WB19 (make sure anhydrous ethanol has been added) to the tube. Mix by vortexing for 15 seconds and proceed to magnetic beads pelleting. Discard the supernatant carefully and avoid pipetting the beads.

5. Repeat step 4 once.

6. Place the microcentrifuge tube on the magnetic stand. Air dry the beads at room temperature until the magnetic beads appear crack and make sure no residual liquid remains in the tube (This step can be performed in a laminar flow cabinet within 10 minutes).

7. Remove the microcentrifuge tube from the magnetic stand and add 100-200 μ l of RNase-free Water to elute DNA/RNA. Mix by vortexing or pipetting up and down for 1 minute and incubate at 65°C for 5 minutes. Gently vortex 2-3 times to resuspend the beads during the incubation.

8. Place the microcentrifuge tube on the magnetic stand for magnetic beads pelleting. Carefully transfer the supernatant into a sterile 1.5 ml microcentrifuge tube and avoid pipetting the beads. Store the eluted DNA/RNA at -70°C.

Notes

- Avoid repeated thawing and freezing samples to ensure high quality DNA/RNA.
- Use RNase-free sterile microcentrifuge tubes and pipette tips to avoid RNA degradation.
- Thoroughly dry beads and the tube before elution to avoid residual ethanol interfering downstream applications.
- For recommendations on automated extraction, please contact Technical Supports from TransGen Biotech.

FOR RESEARCH USE ONLY

